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Nasal Swab Test for selftesting use - Temporary special approval after §11 MPG in Germany (BfArM GZ: 5640-S-179/21)

Please watch the explanatory film before using the self-test on <http://en.wondfo.com.cn/resource/index.html>.

#### INTENDED USE

The 2019-nCoV Antigen Test (Lateral Flow Method) is an immunochromatographic assay for rapid, qualitative detection of novel coronaviruses (2019-nCoV) Nucleocapsid protein extracted from the nasal swab specimen. The test is to be used as an aid in the diagnosis of coronavirus infection disease (COVID-19), which is caused by 2019-nCoV. The test provides preliminary test results. Negative results can not exclude 2019-nCoV infection and they cannot be used as the sole basis for treatment or other management decision. Suitable for ages from 18 years. However, for the users age under 18, it is recommended to perform the test under adult supervision, while the users age over 70 should be aware of the removal of their nasal swab or have nasal swabs assist. For *in vitro* diagnostic use only.

#### SUMMARY

The novel coronaviruses belong to the  $\beta$  genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel corona-virus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

#### PRINCIPLE

The 2019-nCoV Antigen Test (Lateral Flow Method) is an immunochromatographic assay for the rapid, qualitative detection of the antigen of novel coronaviruses (2019-nCoV) obtained from a nasal swab sample. When the extracted specimen is added into the test device, the specimen is absorbed into the device by capillary action, reacts with the 2019-nCoV antibody-dye conjugate and flows across the pre-coated membrane. If the antigen levels of 2019-nCoV are equal to or above the target cutoff, a visible colored band appears in the "test region" (T), indicating that there is a positive result. The absence of this colored band in the test area indicates a negative result. If the test is performed correctly, a line is formed in the "control region" (C) of the test results window. The test provides preliminary test results. Negative test results cannot rule out 2019-nCoV (Corona Disease) infection and cannot be used as the sole basis for treatment or other medical decision.

#### PRECAUTION

- This kit is intended for *in-vitro*-diagnostic use only.
- Unless you are using the test on yourself, it is recommended to wear FFP2 or medical face mask and safety glasses when handling test to the others.
- Proper extraction of nasal swab specimens and accurate performance of the procedure are critical to the result of the test.
- Dispose of all used products after single use with the help of the disposal bag.
- Avoid excessive temperature in the test environment. Test cassette and extraction buffer stored at low temperature must be returned to room temperature before opening in order to avoid condensation of moisture. Do not touch the reaction area of the test strip when opening of the test cassette.
- If the pouch is damaged or the seal is broken, do not use the kit.
- If you have signs of COVID-19 despite a negative result, the test result should be reviewed by a doctor.
- Hold the nasal swab by the style, not the absorbent end.
- This kit is for external use. Don't swallow.
- Avoid getting the buffer solution in contact with eyes or skin.
- Keep out of the reach of children.
- Do not use test kit beyond the expiration date stated on the outer packaging. Always check the expiration date before testing.
- DISPOSAL OF THE DIAGNOSTIC: All components used pose a potential risk of infection. Dispose of them using the disposal bag provided. Stow all the products and materials you have used in it, then close the bag and throw the bag in the trash.

#### STORAGE AND STABILITY

- The products can be stored, transported, stored and used at 2~30°C up to the expiry date printed on the packaging.
- Do not freeze.
- The test cassette must be transported and stored in the sealed pouch.
- The test device should be used within one hour of removing it from the sealed pouch.
- Protect from sunlight, moisture and heat.
- The date of manufacture is printed on the packaging.

#### MATERIALS

		Materials Provided			
Components	REF	W634P0020	W634P0021	W634P0022	W634P0023
Sealed Pouches* (pcs)		1	5	10	20
Drippers (pcs)		1	5	10	20
Extraction Buffer (400 $\mu$ L/tube)		1	5	10	20
Sample Extraction Tube (pcs)		1	5	10	20
Nasal Swab (pcs)		1	5	10	20
Test Tube Rack (pcs)		/	/	1	1
Procedure Card (pcs)		1	1	1	1
Biohazard Bag (pcs)		1	5	10	20
IFU (pcs)		1	1	1	1

Components	REF	W634P0024	W634P0025	W634P0026	W634P0027
Sealed Pouches* (pcs)		1	5	10	20
Pre-installed Extraction Buffer (400 $\mu$ L/tube)		1	5	10	20
Nasal Swab (pcs)		1	5	10	20
Test Tube Rack (pcs)		/	/	1	1
Procedure Card (pcs)		1	1	1	1
Biohazard Bag (pcs)		1	5	10	20
IFU (pcs)		1	1	1	1

**Note:** \*Each sealed pouches containing: 1 Test Cassette and 1 Desiccant Pouch

#### Materials Required but Not Provided

- Timer (Watch, Clock)
- Personal protective equipment (e.g. FFP2 or medical mask, gloves) when handing the contents of this kit for another person.

#### SPECIMEN COLLECTION AND PREPARATION

Bring the test components to room temperature (15 to 30 ° C).

For this test, choose a location where you can run the test as undisturbed as possible for 20 minutes.

Wash or disinfect your hands before performing the testing process.

The test can be performed with a nasal swab sample.

- Take out the nasal swab from the notch.
  - Tilt the head (of the test person) backwards (around 70 degrees).
  - Gently twist the swab, insert the entire absorbent tip of the nasal swab into a nostril for about 1.5 to 2 cm.
  - Perform the first sampling by rubbing the nasal wall firmly with the nasal swab, turning five times circularly against the nasal walls, so that the absorbent surface of the nasal swab is wetted all around.
- Notes:** a) It should not hurt.  
b) Take about 15 seconds to perform this action.
- Slowly remove the nasal swab from the first nostril.
  - Repeat with the same nasal swab in the other nostril.

#### TEST PROCEDURE

##### TEST PROCEDURE I (with Extraction Buffer Tube)


For W634P0020/W634P0021/W634P0022/W634P0023

##### 1. Nasal swab specimen extraction

- Transfer all of the extraction buffer into the sample extraction tube.
- Insert the absorbent nasal swab tip moistened with nasal secretion into the sample extraction tube and twist the nasal swab tip against the bottom and sides of the sample extraction tube ten times to release the secretion from the nasal swab tip. Place the sample extraction tube in the test tube rack and let the swab stand in the sample extraction tube for 1 minute.
- Take the sample extraction tube again in one hand and squeeze it together. With the sample extraction tube squeezed, slowly pull the nasal swab out to extract as much of the liquid as possible from the tip of the nasal swab.
- Place the drippers tip on the sample extraction tube.

##### 2. Test procedure

- Tear open the sealed pouch at the incision, take out the test cassette and place it on a flat surface.
- Turn the sample extraction tube upside down, hold it vertically and add 3, maximum 4 drops processed specimen into the small, round, white sample well. Start the timer. The result window (C T) now turns purple.
- Even if a color line usually appears on control region after approx. 3 minutes, please wait 15 minutes after step (f) and only then read off the test result. This result is valid. After a further 5 minutes, namely more than 20 minutes after step (f), the result is invalid. The results are effective within 15-20 minutes after step (f).

##### TEST PROCEDURE II (with Pre-installed Extraction Buffer)


For W634P0024/W634P0025/W634P0026/W634P0027

##### 1. Nasal swab specimen extraction

- Unscrew the cover of the pre-installed extraction buffer tube. Insert the absorbent nasal swab tip moistened with nasal secretion into the pre-installed extraction buffer tube and twist the nasal swab tip ten times against the bottom and sides of the pre-installed extraction buffer tube to release the secretion from the nasal swab tip.
- Place the nasal swab with the absorbent side in the pre-installed extraction buffer tube. Break off the tip of nasal swab from the first stop (notch) on the end with the swab. The absorbent tip of the nasal swab remains in the pre-installed extraction buffer tube.
- Then unscrew the lid. Place the pre-installed extraction buffer tube in the test tube rack and let the swab stand in the pre-installed extraction buffer tube for 1 minute.

##### 2. Test procedure

- Remove the cap from the lid.
- Tear open the sealed pouch at the incision, take out the test cassette and place it on a flat surface.
- Turn the pre-installed extraction buffer tube upside down, hold it vertically and add 3, maximum 4 drops into the small, round, white sample well. Start the timer. The result window (C T) now turns purple.
- Even if a color line usually appears on control region after approx. 3 minutes, please wait 15 minutes after step (f) and only then read off the test result. This result is valid. After a further 5 minutes, namely more than 20 minutes after step (f), the result is invalid. The results are effective within 15-20 minutes after step (f).

#### WASTE DISPOSAL


All components used pose a potential risk of infection.

- Dispose used products and materials in the biohazard bag.
- Throw the bag in the trash.
- Wash or disinfect your hands

#### RESULT INTERPRETATION

The test contains an internal procedural control. A colored line appearing in the control region(C) is an internal procedural control. It confirms sufficient liquid volume, sufficient moisture penetration of the membrane and correct processing technology.

#### Positive Result

Colored bands appear on both the test line (T) and the control line (C). This indicates a positive result for the 2019-nCoV antigen in the sample. They have most likely become infected with 2019-nCoV. You should immediately go into home isolation (quarantine). The members of your household should also stay away from you immediately. Leave your apartment or house only in medical or other emergencies. Inform your contact persons of the past 14 days about your possible infection. Write down your contacts!

Talk to your doctor about further measures for yourself and your contact persons. Contact the relevant health authority. If necessary, have a PCR test performed to confirm the suspicion of infection with 2019-nCoV. If necessary, inform your employer that you have a positive test result. Observe the quarantine rules! Comply with all the important behavioral and hygiene rules to protect your household members and others from contagion:

- 2m distance (if possible, stay in a separate room)
- Frequent hand washing, other hygiene measures such as disinfection
- Wear medical protective masks or FFP2 masks over the mouth and nose
- Frequent ventilation

If you experience any complaints, seek medical advice immediately. Please note that the legal bases and instructions of the respective responsible government, the district or the city, or the responsible health department apply. There is a possibility of co-infection with other diseases.

#### Negative Result

Only one colored stripe appears on the control line (C). This indicates that the concentration of the 2019-nCoV antigen is zero or below the detection limit of the test. Negative test results do not absolutely rule out an infection. They should not be used as the sole basis for further infection control decisions. Especially in the presence of clinical signs and / or symptoms consistent with COVID-19 (Corona) or for people who have been in contact with the virus. It is recommended that these results be confirmed by a molecular test procedure like PCR test.

#### Invalid Result

No colored band at all appears on control lines (C) and (T), or only a colored band appears on control line (T) after the test has been carried out. The instructions for use may not have been followed correctly or the quality of the test kit may be impaired. It is recommended to repeat the entire process with a new test kit. If the problem persists, stop using the test kits in this package.

#### QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient liquid volume, adequate membrane wicking and correct procedural technique.

#### LIMITATIONS OF PROCEDURE

- The reagent is intended for the detection of 2019-nCoV antigen N protein in human nasal swab samples.
- The accuracy of the test depends on the sampling process. Improper sampling, improper sample storage or repeated freezing and thawing of the sample affect the test result.

- This reagent is a qualitative test. It is not intended to determine the quantitative 2019-nCoV antigen concentration. If you would like to test the quantitative concentration, consult your doctor.
- The test result of the reagent serves only as a clinical reference and should not be used as the sole basis for diagnosis and treatment.
- Due to the limitations of the method of antigen test reagents, it is recommended to use nucleic acid detection or viral culture to verify and confirm negative test results.
- Positive test results do not rule out the possibility of co-infections with other pathogens. A negative result of the reagent can be caused by:

- improper sampling, improper transmission or handling of samples, too low virus titre in the sample;
- The amount of 2019-nCoV antigen is below the detection limit of the test.
- Variations of viral genes can cause changes in antigen determinants.

#### PERFORMANCE CHARACTERISTICS

##### A. Sensitivity and Specificity

356 clinical case samples (including symptomatic and asymptomatic cases) which include 138 confirmed as COVID-19 positive and 218 confirmed as COVID-19 negative by PCR assay, were obtained for testing, and then compared the test results between Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) and the PCR results. The results are shown below.

Reagent	PCR		Total
	Positive	Negative	
Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)	135	2	137
	3	216	219
Total	138	218	356

Sensitivity: 97.83% (95%CI: 93.78%~99.55%)  
Specificity: 99.08% (95%CI: 96.73%~99.89%)  
Total agreement: 98.60% (95%CI: 96.75%~99.54%)

##### B. Cross-reactivity

Cross-reactivity of the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) was evaluated using specimens containing the antigens listed below. The results showed no cross-reactivity with the following:

Common coronavirus antigen (NL63)	1.20*10 <sup>6</sup> TCID <sub>50</sub> /mL	Adenovirus-4 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Common coronavirus antigen (229E)	1.14*10 <sup>6</sup> TCID <sub>50</sub> /mL	Adenovirus-5 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Common coronavirus antigen (OC43)	1.13*10 <sup>6</sup> TCID <sub>50</sub> /mL	Adenovirus-7 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Coronavirus (MERS) antigen	≥10 <sup>6</sup> Copies/mL	Adenovirus-55 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Influenza A H1N1 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL	Enterovirus A/B/C/D antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Influenza A H3N2 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL	EB virus antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Influenza B Yamagata antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL	Measles virus antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Influenza B Victoria antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL	Human Cytomegalovirus antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Respiratory syncytial virus A antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL	Rotavirus antigen	2.76*10 <sup>6</sup> TCID <sub>50</sub> /mL
Respiratory syncytial virus B antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL	Norovirus antigen	3.54*10 <sup>6</sup> TCID <sub>50</sub> /mL
Rhinovirus-A antigen	5.21*10 <sup>6</sup> TCID <sub>50</sub> /mL	Mumps virus antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Rhinovirus-B antigen	2.74*10 <sup>6</sup> TCID <sub>50</sub> /mL	Varicella-zoster virus positive sample	≥10 <sup>6</sup> TCID <sub>50</sub> /mL
Adenovirus-1 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL	Mycoplasma pneumoniae antigen	≥10 <sup>6</sup> Copies/mL
Adenovirus-2 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL	Staphylococcus aureus	3.52*10 <sup>6</sup> CFU/mL
Adenovirus-3 antigen	≥10 <sup>6</sup> TCID <sub>50</sub> /mL		

##### C. Hook effect

Within the titer range of clinically positive samples of 2019-nCoV antigens, there is no hook effect in the test results of this product

##### D. Interference

Can any medication or medical conditions affect the results?

The test result of Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) do not be interfered with the following substance:

Type	Substance	Concentration
Allergic symptoms	Histamine Dihydrochloride	0.25 mg/mL
	Interferon alpha	1.25 mL/L
	Zanamivir	25 mg/L
	Ribavirin	375 mg/L
	Oseltamivir	187.5 mg/L
Antiviral drugs	Palamivir	750 mg/L
	Lopenavir	500 mg/L
	Ritonavir	125 mg/L
	Abidor	0.5 g/L
	Levofloxacin	1.25 g/L
Antibiotics	Azithromycin	2.5 g/L
	Ceftriaxone	2.5 g/L
	Meropenem	500 mg/L
Systemic Antibacterial Drugs	Tobramycin	0.25 mL/L

It is not known whether other drugs will affect the test result. If you have any questions, please contact your doctor or pharmacist and read the information leaflet that came with your medication before using the test.

##### Other possible risks of this test:

- Discomfort during the sampling.
- If you swallow extraction liquid, drink some liquid. If you get them in your eyes, rinse them out with water.

##### E. Precision

- Within run precision was determined by testing positive specimens in 10 times. The agreement rate was 100%.
- Between run precision was determined by testing three different specimens including positive and negative in 3 different lots of test devices. The negative agreement rate and the positive agreement rate were 100%.

##### F. Limit of Detection

The LoD of this test is 1.1×10<sup>2</sup> TCID<sub>50</sub>/mL

#### BIBLIOGRAPHY

- [1] Chen H , Wurm T , Britton P , et al. Interaction of the Coronavirus Nucleoprotein with Nucleolar Antigens and the Host Cell[J]. Journal of Virology, 2002, 76(10).

#### INDEX OF SYMBOL

<b>IVD</b>	In Vitro Diagnostic Use		See Instruction for Use		Expiry Date		Store between 2-30°C
$\Sigma$	Tests per Kit		Manufacturing Date		Keep Dry	<b>REF</b>	Catalog #
<b>LOT</b>	Batch Number		Authorized Representative		Keep away from Sunlight		Do not reuse
	Manufacturer						

Guangzhou Wondfo Biotech Co., Ltd.  
No.8 Lizhishan Road, Science City, Luogang District,  
510663, Guangzhou, P.R.China  
Tel: (+86)400-830-8768  
E-mail: sales@wondfo.com.cn  
Website: en.wondfo.com.cn

**EC REP** Qarad BV  
Cipalstraat 3  
2440 Geel  
Belgium

Test for selftesting use - Temporary special approval after §11 MPG in Germany (BfArM GZ: 5640-S-179/21)